

YK-8000C VITAL SIGNS MONITOR

8 Parameter with 12.1 Inch High Resolution Display



YK-8000C

Features

- 12.1 inch high-resolution color TFT LCD display.
- Thinner and lighter in design, convenient and portable for users.
- Real-time S-T segment analysis, pace-maker detection.
- Various interfaces: standard screen, trend screen, oxyCRG screen, NIBP list screen, big front screen.
- 3-level audio/ visual alarms.
- Efficient resistance to interference of defibrillator and HF knife.
- Up to 400 groups NIBP list, 6000 seconds ECG waveform recall, 60 alarm event records recall, 7-day trend chart in storage.
- Built-in rechargeable lithium battery, 2-hour working capability.
- Wire/ wireless central monitoring system applicable to ICU/ CCU/ OR.

Standard Configuration

SpO₂, PR/HR, NIBP, ECG, RESP, 2-TEMP ETCO₂, 2-IBP,



YK-8000C Vital Signs multi-parameter patient monitor is a robust patient monitor, small enough to move around, yet flexible enough to adapt to your particular needs. With a variety of extras, including arrhythmia and ST analysis, EtCO₂ measurement and double IBP, YK-8000C multi-parameter patient monitor is configurable to suit your requirements. Purchase only the advanced features that you need and nothing more. With NIBP recall and extensive data storage capability, YK-8000C multi-parameter patient monitor is particularly well-suited for use in peri-anesthesia and critical care environments. Rich in functions, conveniently adaptable and pleasingly affordable, YK-8000C multi-parameter patient monitor may be just the monitor for your department.

Features & Benefits:

- 8 parameters: ECG, RESP, NIBP, SPO₂, 2-TEMP, PR/HR, 2IBP, ETCO₂
- 12.1"high resolution color TFT LCD display.
- Real time S-T segment analysis, pace-maker detection.
- Indication of Incorrect manipulation and malfunctionanalysis.
- Patient info input management function
- Multi-lead ECG waveforms display in phase.
- Large volume of tabular and graphic trends informationstorage and easy to recall.
- Capture dynamic waveforms.
- Efficient resistance to interference of defibrillator andHF knife.
- UP to 2 hours working capacity of built-in rechargeable battery.
- Unique accessory management cabinet.
- Networking capacity and build in printer

Support by:



Distributor By :